

Kimberly-Clark* **Critical Care Closed Suction Systems**

KIMBERLY-CLARK* BALLARD* Neonatal and Pediatric TRACH CARE*



As your trusted partner, Kimberly-Clark offers a full range of closed suction systems for your neonatal and pediatric patients.

The KIMBERLY-CLARK* BALLARD* TRACH CARE* Closed Suction Systems are designed to safely suction neonatal and pediatric patients on mechanical ventilation, while providing protection for both the patient and the caregiver.

Our Neonatal and Pediatric Trach Care Closed Suction Systems include:

- **“Soft” low durometer catheter to minimize tracheal trauma**
- **Number and color coded graduations for measured depth suctioning, helping to prevent unnecessary tracheal trauma**
- **Tactile sleeve for catheter advancement**
- **Elbow and “Y” configurations for standard and high-frequency ventilation**
- **Available in 5, 6, 7, 8, 10 & 12 French sizes**

From Medical Devices to surgical apparel – everything we do, everything we make is focused on protecting healthcare professionals and patients. This means we work every day to develop and deliver quality products and to provide responsive customer service that you can depend on to meet your protection needs.

Protection. For life.



Kimberly-Clark* **Critical Care**
Closed Suction Systems

KIMBERLY-CLARK* BALLARD* Neonatal and Pediatric TRACH CARE*

KIMBERLY-CLARK* BALLARD* Neonatal TRACH CARE* Component Kits

Stock #	Outer Diameters	Length	Description	Packaging
1968127	6Fr (2mm)	30.5cm/12 inches	Neonatal "Y" WET PAK* (6 Vials) READY CARE* (5 Ea, 8 Fr Oral Catheters) 9 ft. Suction Tubing	20/Case
1988128	8Fr (2.6mm)	30.5cm/12 inches	Neonatal "Y" WET PAK (6 Vials) READY CARE (5 Ea, 10 Fr Oral Catheters) 9 ft. Suction Tubing	20/Case

KIMBERLY-CLARK* BALLARD* Neonatal and Pediatric TRACH CARE* Elbow

Stock #	Outer Diameters	Length	Description	Packaging
202	5Fr (1.66mm)	30.5cm/12 inches	Elbow Endotracheal Tube Adapters (2mm, 2.5mm)	20/Case
206	6Fr (2mm)	30.5cm/12 inches	Elbow Endotracheal Tube Adapters (2.5mm, 3mm, 3.5mm)	20/Case
207	7Fr (2.3 mm)	30.5cm/12 inches	Elbow Endotracheal Tube Adapters (2.5mm, 3mm, 3.5mm, 4mm)	20/Case
208	8Fr (2.6mm)	35.5cm/14 inches	Elbow Endotracheal Tube Adapters (3mm, 3.5mm, 4mm)	20/Case
210	10Fr (3.3mm)	40.5cm/16 inches	Elbow	20/Case




KIMBERLY-CLARK* BALLARD* Neonatal and Pediatric TRACH CARE* "Y"

Stock #	Outer Diameters	Length	Description	Packaging
195	5Fr (1.66mm)	30.5cm/12 inches	"Y" Endotracheal Tube Adapters (2mm, 2.5mm)	20/Case
196	6Fr (2mm)	30.5cm/12 inches	"Y" Endotracheal Tube Adapters (2.5mm, 3mm, 3.5mm)	20/Case
197	7Fr (2.3mm)	30.5cm/12 inches	"Y" Endotracheal Tube Adapters (2.5mm, 3mm, 3.5mm)	20/Case
198	8Fr (2.6mm)	30.5cm/12 inches	"Y" Endotracheal Tube Adapters (3mm, 3.5mm, 4mm)	20/Case
1910	10Fr (3.3mm)	40.5cm/16 inches	"Y" Endotracheal Tube Adapters (4.5mm, 5mm, 5.5mm)	20/Case
1912	12Fr (4mm)	40.5cm/16 inches	"Y" Endotracheal Tube Adapters (5.5mm, 6mm, 6.5mm)	20/Case




KIMBERLY-CLARK* BALLARD* Neonatal and Pediatric MAC CATHETER

Stock #	Outer Diameters	Length	Description	Packaging
1900	5Fr (1.66mm)	30.5cm/12 inches	MAC Catheter with "Y" Endotracheal Tube Adapters (2.5mm, 3mm, 3.5mm)	5/Case
19020			2.0 "Y" Endotracheal Tube Adapters	50/Case
19040			4.0 "Y" Endotracheal Tube Adapters	50/Case



KIMBERLY-CLARK* BALLARD* Neonatal TRACH CARE* Accessories

Stock #	Description	Packaging
11000	Manual Percussor for Neonatal Pulmonary Use	12/Case
116	Unit Dose Saline 0.9%	12/Case 144/Case



Commitment to Excellence

If, for any reason, our products do not meet your expectations, please let us know your comments or suggestions for improvement. Your input will result in a concerted effort on our part to meet your requirements. Our goal is to provide quality products that completely meet your needs time after time.

For more information, please call 1-800-528-5591 in the United States, +801 572 6800 internationally, or visit our web site at www.kchealthcare.com.

* Registered Trademark or Trademark of Kimberly-Clark Worldwide, Inc. Roswell, GA 30076 or its affiliates. © 2003 KCWW. All rights reserved.

Printed on CLASSIC CREST® Cover, Avon Brilliant White, 65 lb. produced by Kimberly-Clark's Neenah Paper.

Z0413 U6792

Protection. For life.

